

INSTRUCTION MANUAL

ELECTRONIC PERSONAL SCALE

Model: OS-04B

ENGLISH

GB

Dear customers! Thank you for purchasing ORION products. Before operation, please read this instruction manual. Save it for future reference.

IMPORTANT SAFETY INSTRUCTIONS

- This appliance should be used only for its intended purposes. Any other use of it is wrong and therefore dangerous. Manufacturer is not liable for damages that are the result of improper or incorrect use of the appliance. Under these conditions, the guarantee will be cancelled.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Do not immerse scale in water.
- Place the scale in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Strong magnetic fields may interfere with the operation of the scale. Avoid operating the scale in such areas.
- Do not leave any objects on top of scale if the product is not to be used during a prolonged period of time. Also remove the battery to prevent it from being depleted or leaking.
- Do not allow children to play with this appliance.
- Do not use the appliance outdoors.
- The appliance is not intended for commercial use.

USEFUL HINTS

- Place the scale on a hard even surface
- Always weigh yourself on the same scale placed on the same floor surface. Uneven floors may affect the readings.
- Always weigh yourself without clothing or footwear, and before meals. The first thing in the morning is the best time.
- Treat your scale with care as it is a precision instrument, do not drop it or jump up or down on it.
- To avoid slipping, make sure your feet and the glass platform of the scale are dry before standing on the scale.
- To avoid possible injury, do not stand on the edge of the platform.

CHARACTERISTICS

- Equipped with 4 high precision "strain-Gauge" sensor
- LCD Display
- Capacity: 150KG
- Division: 0.1kg
- Auto zero
- Auto off
- Low power & Overload indication
- Power: 1 pc Lithium cell (CR2032)
- Net/Gross weight: 1,35kg/1,70kg
- Unit size : 290x290x20mm
- Gift box size: 303x31x303mm
- Colour: black

PREPARATIONS BEFORE USE

1. Take out your scale.
2. Insert the battery.
3. After you insert the battery the display will show **8888** . Then the display will show 0.0. It means that that the scale has been setup correctly.

Note:

- Make sure that you have installed batteries correctly. " + " is up.

SCALE OPERATION

1. Put the scale on a hard surface (avoiding carpets or any soft surface).
2. Switch on the scale with a slight tap on the glass platform near LCD display.
3. Put your feet evenly on the scale without moving and wait whilst the scale computes weight.
4. After showing the result about 8~12 seconds, the scale will turn off automatically.

Note:

- If you step on the platform before "0.0" appears, you will not obtain a reading.
- If a weight less than 5kg is placed on the scale, the scale will not activate. To weigh something less than 5kg, activate the scale first by tapping the front corner (right or left) of the scale with your foot. Min. weight: 2 kg.

Overload indicator

When the scale is overloaded the following will appear on the screen: ERR

Battery Change Indicator

When the batteries need replacing the following will appear on the screen: LO
Please open the back board of the scale, and then take the batteries carefully out of the battery compartment. Dispose of the old batteries in a safe place. Do not place on fire.

Maintenance

- Wipe dry with a clean cloth
- Do not use abrasive washing powders, steel wool or scouring pads
- Store scale in a cool place
- Do not store the scale in an upright position